# Weighted kernel density estimation to quickly reproduce the profile of a diffractometer

This example shows a work-arround for a quick visualization of a diffractorgram (similar to experimental powder diffractograms) from **ImageD11** ".flt" or ".new" columnfile containing peaks information.

It is basically a **probability density function** (pdf) of the  $2\theta$  *position* of the peak, which is weighted by the *peak intensity*.

The smoothing of such gaussian kde is decided by the *bandwidht* value.

<u>Weighted kde (https://nbviewer.jupyter.org/gist/tillahoffmann/f844bce2ec264c1c8cb5)</u>: The original Scipy gaussian kde was modified by <u>Till Hoffmann (http://tillahoffmann.github.io/)</u> to allow for heterogeneous sampling weights.

Loading and visualizing the input data


```
Entrée [3]: # read the peaks
flt = columnfile('sma_261N.flt.new')

# peaks indexed to phase 1
phase1 = flt.copy()
phase1.filter( phase1.labels > -1 )

# unindexed peaks (phase 2 + unindexed phase 1?)
phase2 = flt.copy()
phase2.filter( phase2.labels == -1 )

#plot radial transform for phase 1
plt.plot( phase1.tth_per_grain, phase1.eta_per_grain, 'x')
plt.xlabel( r'$ 2 \theta \, (\degree) $' )
plt.ylabel( r'$ \eta \, (\degree) $' )
plt.title( r'$Diffraction \, angles$' )
```

Out[3]: Text(0.5, 1.0, '\$Diffraction \\, angles\$')




#### Plotting the diffraction profile

```
Entrée [4]: # Probability density function (pdf) of 2theta
    # weighted by the peak intensity and using default 2theta bandwidth
    I_phase1 = phase1.sum_intensity * phase1.Lorentz_per_grain
    pdf = wkde.gaussian_kde( phase1.tth_per_grain, weights = I_phase1)

# Plotting it over 2theta range
    x = np.linspace( min(flt.tth), max(flt.tth), 500 )
    y = pdf(x)
    plt.plot(x, y)
    plt.xlabel( r'$ 2 \theta \, (\degree) $' )
    plt.ylabel( r'$ I $' )
    plt.yticks([])
    plt.title( ' With bandwidth = %.3f'%pdf.factor )
```

Out[4]: Text(0.5, 1.0, ' With bandwidth = 0.247')

#### With bandwidth = 0.247



The profile showed above is highly smoothed and the hkl peaks are merged.

→ A Smaller bandwidth should be used.

### Choosing the right bandwidth of the estimator

The bandwidth can be passed as argument to the gaussian\_kde() object or set afterward using the later set\_badwidth() method. For example, the bandwidth can be reduced by a factor of 100 with respect to its previous value:

```
gaussian kde().set bandwidth( gaussian kde().factor / 100 )
```

```
Entrée [5]:
            pdf phase1 = wkde.gaussian kde( phase1.tth, weights = phase1.sum intensity )
            pdf_phase2 = wkde.gaussian_kde( phase2.tth, weights = phase2.sum_intensity )
            frac_phase1 = np.sum( phase1.sum_intensity ) / np.sum( flt.sum_intensity )
            frac_phase2 = np.sum( phase2.sum_intensity ) / np.sum( flt.sum_intensity )
            from ipywidgets import interact
            bw_range = ( 0.001, pdf_phase1.factor/3, 0.001)
            @interact( bandwidth = bw_range)
            def plot_pdf(bandwidth):
                pdf phase1.set bandwidth(bandwidth)
                pdf phase2.set bandwidth(bandwidth)
                y_phase1 = pdf_phase1(x)
                y_phase2 = pdf_phase2(x)
                plt.plot( x, frac_phase1 * y_phase1, label = r'$Phase \, 1$')
                plt.plot( x, frac_phase2 * y_phase2, label = r'$Phase \, 2$' )
                plt.legend(loc='best')
                plt.xlabel( r'$ 2 \theta \, (\degree) $' )
                plt.ylabel( r'$ I $' )
                plt.yticks([])
                plt.title( r'$ 3DXRD \, diffractogram $' )
```

bandwidth 0.01

## $3DXRD\, diffractogram$

